

SIEMENS

MULTIMOBIL 5C

SP

Pre-Installation

© Siemens AG 2004

The reproduction, transmission or use of this document or its contents is not permitted without express written authority. Offenders will be liable for damages. All rights, including rights created by patent grant or registration of a utility model or design, are reserved.

English

Doc. Gen. Date: 10.05

Print No.: SPR2-125.811.01.02.02

Replaces: SPR2-125.811.01.01.02

Multimobil 5C

Med



Pre-installation

Version 5.0

Copyright © SIEMENS LTD., Med INDIA

1	PRE-INSTALLATION	3
1.1	TECHNICAL SPECIFICATIONS	3
1.2	INFORMATION ON ROOM PLANNING	7
1.3	PACKAGING	7
1.4	UNLOADING	8
1.5	UNPACKING THE MM5C	8
1.6	UNPACKING THE MONITOR TROLLEY	10
1.7	INSTALLATION OF THE MONITORS	11

1 Pre-installation

1.1 Technical Specifications

Sr.No	Parameter	Value
1.	Mains Voltage	110V/190 -240V~
2.	Line frequency	50 / 60Hz
3.	Line impedance	1.5 Ω (max) @ 190-240 V 0.36 Ω (max) @ 110V
4.	Length of power cable	5 m
5.	Power Output – Radiography Nominal Electric power at 100kV and 100mSec.	2.5 kW (100 kV, 25 mA) @ 190-240V 1.6KW (100 kV, 16 mA) @110V
6.	Wave Shape	Multipulse DC - Ripple 5kV max
7.	KVp Range	40 – 100kV in 20 steps
8.	KV Accuracy	$\leq \pm 5 \%$
9.	mA-range	13 – 63 mA @ 190-240V input 10 – 40 mA @ 110V input
10.	mAs Range	0.32 – 200mAs at 40kV 0.32 – 160mAs at 42 - 50kV 0.32 – 125mAs at 52 - 63kV 0.32 – 100mAs at 66 - 77kV 0.32 – 80 mAs at 81 - 100kV
11.	mAs Accuracy	$\leq 10 \%$ + 0.2 mAs for mAs ≤ 20 mAs $\leq 5 \%$ + 0.2 mAs for mAs > 20 mAs
12.	Exposure Time	20mSec - 5 Sec in 24 steps @190-240V 20mSec – 8 Sec in 26 steps @ 110V
		For details of kV, mA and exposure time combinations as per EN 60601-2-7 refer to exposure chart on page 6-14.
13.	Fluoroscopy voltage	40-110 kV
14.	Fluoroscopy Modes	ADR(Automatic dosage regulation): Anti-isowatt & Contrast.
15.	Fluoroscopy mA range	Anti-isowatt 0.5 to 5 mA Contrast 0.5 to 7mA
16.	Fluoroscopy Power	550W continuous
17.	Timer	Fluoroscopy time display with timer

Sr.No	Parameter	Value
		reset and 5 minute buzzer indication
18.	Thermal Protection	Auto cutoff of fluoroscopy and radiography in case of tube head temperature exceeding 60°(+5°) Celsius
19.	Automatic termination	Buzzer alarm is set and fluoroscopy terminates after 10 minutes of Continuous uninterrupted fluoroscopy.
20.	Power Input: Continuous Instantaneous	1.5kVA 3.5kVA
21.	Mains Isolation External: done by the user	Power supply cord shall be plugged where both poles (L & N) are isolated simultaneously.
22.	Type and degree of protection against electric shock (EN 60601 – 1)	Class – I, Type B equipment
23.	X - ray Tube	Stationary Anode DF-151R
24.	Focal Spot – nominal value	0.5 for Fluoroscopy 1.5 for radiography
25.	Inherent tube Filtration	0.8 mm Al at 50kV
26.	Nominal kV	110 kV
27.	Target Angle	16°
28.	Collimator	Motorised, multileaf Iris Collimator
29.	Total filtration	5.3 mm Al at 100kV
30.	Exposure release Switch. Flouroscopy switch.	Hand held - 1 Step. Foot switch – 1 step and / or Hand held switch – 1 step.
31.	Exposure Rate for Radiography	Pulse to pause ratio 1:30; corresponds to a cool down period of 4 minutes at maximum output.
32.	Max. X-ray coverage	430X430mm at SID 909 mm
33.	Mode of Operation	Continuous operation
34.	Focus-II distance	860mm for 9" system 900mm for 6" system
35.	Max. cassette size at 900 cm SID	10"X12"

Sr.No	Parameter	Value
36.	Grid	Circular grid placed at the input of the II tube , 40 lp/cm, 8:1
37.	Image Intensifier	9" / 6"
38.	Image Intensifier tube resolution	Nominal 52lp/cm and zoom 68 lp/cm for 9" system Nominal 66 lp/cm and zoom 77 lp/cm
39.	Nominal entrance field	230mm for 9" system 150mm for 6" system
40.	II output image size diameter	25±0.5mm for 9" system 20±0.5mm for 6" system
41.	Zoom (field size)	230mm/160mm for 9" system 150 mm/100 mm for 6" system
42.	TV camera Input Voltage	AC100~220V±10%
43.	Pick Up device Pick UP area	CCD 6.4X4.8mm
44.	Effective Pixels	752(H)X582(V)(CCIR) 768(H)X494(V)(EIA)
45.	TV System	625lines/50Hz(CCIR) 525lines/60Hz(EIA)
46.	No. of monitors	2 monitors one each for Live/Last Image Hold(LIH) and stored memory display
47.	Monitor size	17"monochrome with 9" system and 15" monochrome with 6" system
48.	Monitor screen	High brightness, high resolution with independent brightness and contrast control
49.	TV standard	CCIR 625 lines at 50 Hz
50.	Aspect ratio	4:3
51.	Image rotation	Endless clockwise and anti-clockwise with complete 400°rotation
52.	Recursive filter	Steps of OFF,2,4,8,16 to minimize dynamic blur
53.	Image reversal	Left/right and top/bottom
54.	II Chain Resolution	Nominal 1.2 lp/cm Zoom 1.6 lp/cm for 9" system Nominal 1.4 lp/cm Zoom 2.0 lp/cm for 6" system

Sr.No	Parameter	Value
55.	C-arm orbital movement	125°(+90°to -35°)
56.	Angulation	±190°
57.	Horizontal travel	200mm
58.	Swivel range	±12.5°
59.	Vertical travel	400 mm
60.	Mobility	Directional lock on rear wheels enabling excellent maneuverability
61.	Max. floor incline for transport	5°
62.	Mechanical Dimensions (lxbxh)	2000 x 800 x1720 mm (min) 2070x1050 x2120 mm (max)
63.	Weight	
	Without packing	280kg
	With packing	680kg
64.	Environmental conditions Transport and Storage Temp range.	-15° C to 45° C
	Relative Humidity range	Upto 95% (No Condensation)
	Atmospheric pressure	76 kPa to106 kPa
64.	Operating Temperature	10° C to 40° C
	Relative Humidity range	Upto 85% Max non condensing
	Atmospheric pressure	76 kPa to106 kPa
65.	Conformance to Standards	EN 60601-1:1990 EN 60601-2-7:1998 EN 60601-1-3:1994 EN 60601-1-2:2001 Compliance to AERB Type Approval Certification Compliance to BIS test certification.

1.2 Information on Room Planning

No specific conditions required. Following are the recommended room conditions for proper functioning of the unit.

Ambient Temperature (Room temperature): + 40 °C Max

Air Humidity: up to 75%, Non-condensing.

During Transport and storage of the product the ambient temperature should not be below -15°C or rise above +45°C. Storage is permissible in rooms with minimum amount of dust and Humidity in range of 30% to 95% RH provided no condensation occurs.

1.3 Packaging

The Multimobil 5C unit is delivered in two wooden boxes.

One box consists of the Monitor trolley. This is a wooden box with the trolley mounted on the wooden base and is fixed using clamping brackets.

The second box contains the MM5C unit along with monitors and accessories .The C-arm unit is clamped using brackets on the wooden base. The two monitors are packed in the corrugated box and placed on the wooden platform.

One box consists of standard accessories i.e. sterile cloth, Clips and user manual.



MM5C with
monitors and
accessories



Monitor trolley
for MM5C

1.4 Unloading

The unloading of the unit from the truck can be carried out using any of the following ways depending on the site conditions:

- 1) Using a forklift.
- 2) Using an overhead travelling crane.
- 3) Using a chain pulley block and tripod.
- 4) Using conventional ways i.e. using wedges, wooden planks, pipes or rods and labour etc.

Note:

1. All boxes are provided with transport base.
2. Do not roll the boxes.
3. Do not drop the boxes from a height of more than 20cms.
4. Observe the stickers **(This side up)** and **(handle with care)** on the Packing for proper handling.

Once the unit is removed from the transport, place it on the ground on the base.


1.5 Unpacking the MM5C

The complete system must be inspected for physical damage and / or material short-fall and factory must be informed accordingly.

The Multimobil 5C consists of fragile components like camera, Image intensifier, X-ray Tube etc. Sufficient care must be taken while unloading / unpacking the unit. **Refer the unpacking instructions provided on the packing.**

Observe the following steps:

- 1) Remove the restraining strips by cutting it. Use a cutting plier for this purpose.
- 2) Start removing the wooden planks. For this purpose the service engineer should use a claw hammer or nail pulling plier and preferably a crow bar.
- 3) First of all remove the top planks. Upon opening this plank you will be able to see the unit inside. Ensure that the unit has not received any transit damages.
- 4) Now, remove the front and rear wall of the box one by one. The front wall is provided with a ramp. Hence, place this plank carefully for pushing the unit out of the platform.
- 5) Now you will see three wooden blocks fixed across the unit for supporting the wall and also to support the image intensifier and single tank. These are supporting beams.
- 6) Loosen these blocks from one side and remove them gently. Now, remove the side walls.
- 7) Care should be taken to see that the wooden planks or the nails do not damage the unit.
- 8) Now, you can easily remove the ESD bags.
- 9) After removing the ESD bags, remove the monitor boxes and accessories and take them to the OT.
- 10) Remove the packing strips and thermocol packing from the control, II and single tank. Also, remove the wooden lever fixed at the base by cutting the strap.
- 11) Next step is to remove the clamping brackets. First of all, remove the bracket used for clamping the single tank and Image intensifier. Secondly, remove the bracket used to clamp C-arm. Third, remove the brackets on the rear wheels. Next remove the brackets used to clamp the base.
- 12) Now ensure that the rear wheels are in locked condition. Also the C-Arm should be in locked condition.
- 13) Remove the bolts of the clamping brackets at the rear side on which the housing is resting. Now, using the lever, gently lift the base along the bracket so that the rubber block can be removed. Once the rubber blocks are removed, remove the clamping bracket.
- 14) Now, remove the two brackets on the front end of the base. Lift from sides using wooden lever. Remove the rubber buffer and remove the bracket.
- 15) Now the unit will be resting on the wooden platform. Release the brakes and push the unit down to the ground through the ramp.



DO NOT HOLD THE UNIT AT POINTS

MARKED  WHILE LIFTING

Observe the following steps for transporting the unit to the OT.

Unlock the C-arm and turn it to the position shown in the figure and lock it. This is the suitable position for moving the unit from the platform.

Now, release the brake of the C-Arm and swing the C-Arm back to the normal position as shown in the figure for transportation. This is also the parking position.

Roll all the cables and hang it on the cable holder provided on the housing.

Unbrake the wheels and move the unit to the OT. Use the steering for easy manoeuvrability.

1.6 Unpacking the monitor trolley

Observe the following steps:

- 1) Remove the restraining strips by cutting it .Use a cutting plier for this purpose.
- 2) Remove the top plank.
- 3) Remove the front and rear wall of the box.
- 4) Remove the supporting beams fixed across the sidewalls.
- 5) Remove the side planks. The trolley is mounted on the clamping brackets.
- 6) Remove the ESD bag.
- 7) Now, using a spanner of 16-17, remove the clamping brackets.
- 8) Lift the trolley gently by holding at the base and supporting it from top. Place it on the floor.
- 9) Release the brakes and the trolley can be shifted to the OT.

1.7 Installation Of The Monitors

- 1) Open the corrugated box and remove the monitors.
- 2) Open the accessories box and remove the packet of hardware.
- 3) Place the monitors on the top tray of the monitor trolley and fasten them using socket HD screw M4x20, spring washer and plain washers.
- 4) Make use of the connector provided on the monitor trolley for providing power supply to the monitors.
- 5) Ensure that the radiation indication lamp is available on the Monitor Trolley.